



SOFT NEXT

The new bandage generation











16 16 16 16 16











FACTS A	ND BE	NEFITS
---------	-------	--------

ADHESIVE-FREE BANDAGE	Does not stick to skin, hair nor wounds No glue residue
ONLY ADHERES TO ITSELF	Can be used directly on wet and dirty skin. Can be used on top of creams and ointments
LATEX FREE	Skin and allergy-friendly
SOFT AND ELASTIC	Shapes itself to the joints Fits under gloves
ABSORBENT	Ideal for bleeding wounds
STOPS BLEEDING	Place Soft NEXT over the wound and pinch the bandage until the desired pressure is reached. Use several layers if the bleeding is severe
CAN BE USED IN WATER	Does not fall off when wet
POROUS	Does not leave wrinkled skin beneath the plaster
SIZES: 6 cm x 100 cm 6 cm x 450 cm 3 cm x 450 cm	Different sizes cater to different needs
	The blue version is highly visible and well

NEUTRAL OR BLUE

The blue version is highly visible and well suited for the food industry





Easy to use • Only adheres to itself • Leaves no glue residue • Soft and elastic • Shapes itself to the joints • Can be used in water

Snøgg AS • Postboks 5444 Strai, NO-4671 Kristiansand Tlf: +47 38 03 90 60 • Faks: +47 38 03 90 61 E-mail: snogg@snogg.no **Distributor:**

snogg.no | snogg.com